Amplatz Type Renal Dilator/Sheath Set

en	Instructions for Use	2
es	Instrucciones de uso	7
fr	Mode d'emploi	12
de	Gebrauchsanweisung	17
it	Istruzioni per l'uso	22
nl	Instructies voor gebruik	27
pt-EU	Instruções de Utilização	32

51649588-01

2024-12

TABLE OF CONTENTS

REUSE WARNING3

DEVICE DESCRIPTION3

 Contents.....3

 Operating Principle3

 User Information.....3

INTENDED USE/INDICATIONS FOR USE3

 Clinical Benefit Statement.....4

CONTRAINDICATIONS.....4

WARNINGS.....4

PRECAUTIONS.....4

ADVERSE EVENTS4

HOW SUPPLIED.....4

 Device Details.....4

 Handling and Storage.....4

OPERATIONAL INSTRUCTIONS5

 Procedure.....5

 Disposal5

 Post-Procedure5

PATIENT COUNSELING INFORMATION5

WARRANTY.....6

SYMBOL DEFINITIONS.....6

Amplatz Type Renal Dilator/Sheath Set

Rx ONLY

Caution: Federal Law restricts this device to sale by or on the order of a physician.

REUSE WARNING

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

DEVICE DESCRIPTION

The Amplatz Type Renal Dilator/ Sheath Set is used to progressively dilate the nephrostomy tract during percutaneous nephrostomy procedures and to maintain the tract opening once established to allow passage of surgical instrumentation. The Set includes 12 Amplatz Type Renal Dilators 8F through 30F (even sizes) of 35 cm length, 4 Amplatz Type Renal Sheaths 24F through 30F (even sizes) of 17 cm length, and an 8/10 Dilator/Sheath Set that includes one 8F stylet (70 cm long) with an 8F ID/10F OD radiopaque polytetrafluoroethylene (PTFE) sheath (30 cm long). The Amplatz Type Dilators (35 cm long) have a gradual tapered tip designed for atraumatic placement and manoeuvrability. The Amplatz Type Renal Sheaths are made of firm, radiopaque PTFE with tapered tips for atraumatic placement and ease of advancement over the dilators.

Contents

- (12) Dilators, 8F (2.7 mm) - 30F (10.1 mm) (Even sizes), 35 cm Length
- (4) Sheaths, 24F (8.1 mm) - 30F (10.1 mm) (Even sizes)
- (1) Stylet, 8F (2.7 mm) O.D., 70 cm Length
- (1) Sheath, 8F (2.7 mm) I.D., 10F (3.4 mm) O.D. Radiopaque PTFE, 30 cm Length

Operating Principle

The Amplatz Type Renal Dilator/Sheath Set is used to progressively dilate and maintain the nephrostomy tract opening once initial access has been established.

User Information

This device is restricted to use by or under the supervision of physicians trained in nephrostomy procedures.

INTENDED USE/INDICATIONS FOR USE

The Amplatz Type Renal Dilator/Sheath Set is recommended for dilatation of the nephrostomy tract and sheath placement.

Clinical Benefit Statement

The clinical benefit of the Amplatz Type Dilator/ Sheath Set is to progressively dilate the nephrostomy tract for sheath placement, in order to provide a percutaneous access lumen through which surgical instrumentation can be passed for therapeutic or diagnostic procedures.

CONTRAINDICATIONS

This product is contraindicated in the presence of conditions which create unacceptable risk during catheterization.

WARNINGS

None known.

PRECAUTIONS

A thorough understanding of the technical principles, clinical applications, and risks associated with nephrostomy dilatation is necessary before using this product.

Do not advance the Amplatz Type Dilator/ Sheath if significant resistance is encountered. Forceful insertion may result in patient injury.

The recommendations given are meant to serve only as a basic guide to the utilization of this set. The performance of percutaneous nephrostomy should not be undertaken without comprehensive knowledge of the indications, techniques and risks of the procedure.

ADVERSE EVENTS

The adverse events which may result from a dilatation procedure include:

- Avulsion
- Bleeding
 - o Hematoma
 - o Hematuria
 - o Hemorrhage
- Edema
- Extravasation
- Fistula
- Infection
- Laceration
- Obstruction
- Pain
- Renal Insufficiency/ Failure
- Tissue perforation
- Tissue trauma
- Unretrieved Device Fragments (UDF)

Adverse events may require additional medical or surgical intervention.

HOW SUPPLIED

Device Details

The Amplatz Type Renal Dilator/Sheath Set is supplied sterile for single use using an ethylene oxide process. The devices are arranged in a tray and the tray is placed into a pouch and sealed which acts as the sterile barrier. Do not use if the sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

- Do not use if package is damaged or unintentionally opened before use.
- Do not use if labeling is incomplete or illegible.
- Do not use past the "Use By" date displayed on the product label.

Handling and Storage

This product has no special handling or storage requirements.

OPERATIONAL INSTRUCTIONS

Procedure

1. Establishment of a nephrostomy tract and renal dilatation should always be performed under appropriate imaging guidance (fluoroscopy).
2. Establish percutaneous access and place an 0.038 in (0.97 mm) diameter guidewire into a secure location within the intra-renal system. If feasible, pass the guidewire down the ureter and into the bladder.
3. Progressively dilate the tract over the guidewire with the 8F through 12F (2.7 mm - 4.1 mm) dilators while holding the guidewire in place.
4. Pass the 8F (2.7 mm) catheter/stylet and sheath over the guidewire into the tract.
5. Remove the 8F (2.7 mm) catheter/stylet while maintaining the position of the guidewire and sheath.
6. Place additional 0.038 in (0.97 mm) diameter guidewire(s) through the sheath into a secure location in the intra-renal system.
7. Remove the sheath maintaining the position of the guidewires. Lay the safety guidewire(s) aside and affix it to the drape.
8. Pass the 8F (2.7 mm) catheter/stylet over the working guidewire.
9. The tract may now be progressively dilated from 14F through 30F (4.7 mm - 10.1 mm) as needed by passing the dilators over the 8F (2.7 mm) catheter/stylet.
10. When dilating the tract, hold the 8F (2.7 mm) catheter/stylet in position and pass the dilators into the tract using a twisting motion.
11. Once the tract has been dilated to the desired size, pass the appropriately sized sheath over the dilator and position as desired.
12. Remove the 8F (2.7 mm) catheter/stylet, maintaining the position of the guidewire, dilator and sheath.
13. Remove the dilator from the sheath, maintaining the position of the sheath and guidewire.
14. Appropriate instrumentation can now be passed through the sheath.

Disposal

To minimize the risk of infection or microbial hazards after use, dispose of the device and packaging as follows:

After use, the device and packaging may contain biohazardous substances. Any device and packaging that came into contact with biohazardous substances should be treated and disposed of as biohazardous waste or be treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with a biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

Post-Procedure

Any serious incident that occurs in relation to this device should be reported to Boston Scientific and the relevant local regulatory authority.

PATIENT COUNSELING INFORMATION

Brief the patient on any relevant post-procedure instructions, contraindications, warnings, precautions, and/or adverse events found within this Instructions for Use pertaining to the patient. Inform the patient if any serious incident occurs in relation to this device, it should be reported to Boston Scientific, and to the relevant local regulatory authority.

WARRANTY

For device warranty information, visit www.bostonscientific.com/warranty

EU Importer: Boston Scientific International B.V., Vestastraat 6, 6468 EX
Kerkrade, The Netherlands

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary.

The following symbol is used for Contents

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or its affiliates. All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
Boston Scientific Corporation
© 2024 Boston Scientific Corporation