RescueNet™

Retrieval Device

en	Instructions for Use	2
es	Instrucciones de uso	7
fr	Mode d'emploi	12
de	Gebrauchsanweisung	17
it	Istruzioni per l'uso	22
nl	Instructies voor gebruik	27
pt-EU	Instruções de Utilização	32

52084303-01

2025-04

TABLE OF CONTENTS

REUSE WARNING3

DEVICE DESCRIPTION3

 Figure 1. RescueNet Retrieval Device3

 Contents3

 Operating Principle3

 Materials3

 User Information3

INTENDED USE/INDICATIONS FOR USE4

 Clinical Benefit Statement4

CONTRAINDICATIONS4

WARNINGS4

PRECAUTIONS4

ADVERSE EVENTS4

HOW SUPPLIED4

 Device Details4

 Handling and Storage4

OPERATIONAL INSTRUCTIONS5

 Disposal5

 Post-Procedure5

PATIENT COUNSELING INFORMATION5

WARRANTY6

SYMBOL DEFINITIONS6

 Table 1. Symbol Definitions6

RescueNet™

Retrieval Device

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

DEVICE DESCRIPTION

The device consists of a flexible wire cable and snare loop with net which can be extended and retracted from the flexible outer sheath using a three ring handle. The device is for use with an endoscope with a working channel that is 2.8 mm or greater.


Figure 1. RescueNet Retrieval Device

Contents

Each packaged pouch contains the following:

- (1) RescueNet Retrieval device

Operating Principle

The RescueNet device is a sterile device comprised of a flexible wire cable and snare loop with net, which can be extended and retracted from the flexible outer sheath using a three ring, proximal handle. Objects are retrieved by advancing the finger rings of the handle to open the retracted net. The net is then endoscopically manipulated over the foreign object, food bolus, tissue fragment or excised tissue, in which the finger rings of the handle are retracted to close the net. The RescueNet device and endoscope are removed from the patient together once the article has been retrieved.

Materials

Contains cobalt: CAS No. 7440-48-4; EC No. 231-158-0. Defined as a 1B carcinogen and reproductive toxicant according to the European Commission in a concentration above 0.1 % weight by weight.

Note: This device is made with a metal alloy which contains cobalt. Current scientific evidence supports that metal alloys containing cobalt used in medical devices do not cause an increased risk of cancer or adverse reproductive effects.

User Information

Read the entire Instructions for Use thoroughly before using RescueNet. This device should only be used by or under the supervision of physicians thoroughly trained in endoscopic retrieval techniques. A thorough understanding of the technical principles, clinical applications, and risks associated with endoscopy procedures is necessary before using the device.

INTENDED USE/INDICATIONS FOR USE

This device is intended for use in endoscopic retrieval of foreign objects, food bolus, tissue fragments and excised tissue such as polyps.

Clinical Benefit Statement

Endoscopic retrieval of foreign objects, food bolus, tissue fragments and excised tissue such as polyps.

CONTRAINDICATIONS

Contraindications include those specific to the primary endoscopic procedure to be performed in gaining access to the desired site.

WARNINGS

- Do not use this device for any purpose other than the stated intended use.
- Avoid blindly passing the device past any foreign object, particularly if the entire lumen is blocked to avoid patient injuries.
- Use care to avoid inadvertently grasping tissue or organs not intended for retrieval to avoid patient injuries such as bleeding, tissue damage, etc.
- Keep gentle traction on the device during retrieval so that the retrieved object does not become loosened, dislodged or aspirated into the trachea causing injuries such as obstruction, respiratory distress, and/or inflammation.
- Do not use this device for the removal of sharp objects as this can cause injuries such as perforation, bleeding, and/or tissue damage. Sharp objects may also damage the net and could cause the retrieved object to dislodge inside the patient.

PRECAUTIONS

- After extraction, if the retrieved material does not exit the net upon deployment, the material should be removed from the net by rinsing the net in a basin of water; do not "pick" retrieved material from the net with a finger or tool as this may damage the net.

ADVERSE EVENTS

- Aspiration
- Bleeding
- Infection (including pneumonia)
- Inflammation/Edema
- Obstruction
- Pain
- Perforation
- Respiratory Distress/Insufficiency
- Tissue Damage

HOW SUPPLIED

Contents supplied STERILE using an ethylene oxide (EO) process. Before using, carefully examine the unit to verify that neither the contents nor the sterile package have been damaged in shipment. Do not use if damaged. Immediately return damaged product to Boston Scientific Corporation.

Device Details

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

Handling and Storage

This product has no special handling or storage requirements.

Check device label for expiration date.

OPERATIONAL INSTRUCTIONS

1. Prior to advancing the device into the endoscope, uncoil the device and drape in a “U” shaped configuration. The device can be damaged if it is coiled while the handle is operated.
2. Holding the proximal end of the device in one hand and the distal sheath in the opposite hand, slowly retract the finger rings until the net is completely withdrawn into the sheath. Confirm that there are no significant deformations or abnormalities.
3. Once the foreign object, food bolus, tissue sample or excised polyp for retrieval has been identified and with the net fully retracted within the catheter, advance the catheter into the accessory channel using short, 1.0 cm – 1.5 cm strokes until the distal end of the sheath is endoscopically visualized.
4. Advance the distal end of the sheath slightly past the foreign object, food bolus, tissue sample or excised polyp, however, avoid blindly passing the device past any foreign object, food bolus, tissue sample or excised polyp if the entire lumen is blocked. Open the net by advancing the finger rings forward. Confirm that the net is fully open via endoscopic observation.
5. Endoscopically manipulate the net over the foreign object, food bolus, tissue sample or excised polyp.
6. Retrieve the foreign object, food bolus, tissue sample of excised polyp by slowly retracting the finger rings on the handle. Gentle traction should be applied on the handle to keep the net closed.
7. Retract the catheter retrograde toward the endoscope so that the net is endoscopically observed yet not obscuring the endoscopic view. Endoscopic observation is necessary during endoscope and device removal so as not to lose sight of the foreign object, food bolus, tissue sample or excised polyp during removal.

Precaution: Remove the endoscope and net simultaneously. Do not attempt to pull the net and captured object through the working channel.

8. Once the endoscope and the net have been removed, the retrieved foreign object, food bolus, tissue sample or excised polyp can be removed from the net by advancing the finger rings forward to open the net. Specimens should be prepared according to standard technique for histologic evaluation.

Disposal

To minimize the risk of infection or microbial hazards after use, dispose device and packaging as follows:

After use, device may contain biohazardous substances. The device and packaging should be treated and disposed of as biohazardous waste or have them treated and disposed of in accordance with any applicable hospital, administrative, and/or local government regulations. Use of a biohazardous container with biological hazard symbol is recommended. Untreated biohazardous waste should not be disposed of in the municipal waste system.

Post-Procedure

Any serious incident that occurs in relation to this device should be reported to the manufacturer, and to the relevant local regulatory authority.

PATIENT COUNSELING INFORMATION

Brief the patient on any relevant post-procedure instructions, contraindications, warnings, precautions and/or adverse events found within this Instructions for Use (IFU), pertaining to the patient.

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

The following is a trademark of Boston Scientific Corporation or its affiliates: RescueNet.

All other trademarks are the property of their respective owners.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary.

Table 1. Symbol Definitions

	Contents
	Minimum Required Working Channel

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or its affiliates. All rights reserved.
A seguinte marca é uma marca comercial da Boston Scientific Corporation ou das respetivas afiliadas: RescueNet.
Boston Scientific Limited
© 2025 Boston Scientific Corporation