# SpyGlass™ DS Digital Controller

## User's Manual

---

|                              |     |                                  |     |
|------------------------------|-----|----------------------------------|-----|
| User's Manual .....          | 2   | Käyttöopas .....                 | 275 |
| Brugervejledning .....       | 26  | Manual de utilizar .....         | 299 |
| Εγχειρίδιο χειρισμού .....   | 51  | Používateľská príručka .....     | 324 |
| Manual do Usuário .....      | 76  | Ръководство за потребителя ..... | 349 |
| Bruksanvisning .....         | 101 | Korisnički priručnik .....       | 374 |
| Felhasználói kézikönyv ..... | 125 | Kasutusjuhend .....              | 399 |
| Uživatelská příručka .....   | 150 | Lietotāja rokasgrāmata .....     | 423 |
| Podręcznik użytkownika ..... | 175 | Naudotojo vadovas .....          | 447 |
| Brukerhåndbok .....          | 200 | Pedoman Pengguna .....           | 472 |
| 사용 설명서 .....                 | 225 | Uporabniški priročnik .....      | 497 |
| Kullanıcı Kılavuzu .....     | 250 |                                  |     |

|                                                                               |           |
|-------------------------------------------------------------------------------|-----------|
| <b>TABLE OF CONTENTS</b>                                                      |           |
| <b>INTRODUCTION</b>                                                           | <b>3</b>  |
| Where to Get Help                                                             | 3         |
| <b>DEVICE DESCRIPTION</b>                                                     | <b>3</b>  |
| Contents                                                                      | 3         |
| User Information                                                              | 4         |
| <b>INTENDED USE/INDICATIONS FOR USE</b>                                       | <b>4</b>  |
| Clinical Benefit Statement                                                    | 4         |
| <b>CONTRAINDICATIONS</b>                                                      | <b>4</b>  |
| <b>WARNINGS</b>                                                               | <b>5</b>  |
| <b>PRECAUTIONS</b>                                                            | <b>6</b>  |
| <b>ADVERSE EVENTS</b>                                                         | <b>6</b>  |
| <b>CONFORMANCE TO STANDARDS</b>                                               | <b>6</b>  |
| Essential Performance Statement                                               | 6         |
| <b>HOW SUPPLIED</b>                                                           | <b>6</b>  |
| Handling and Storage                                                          | 6         |
| Environmental Limits During Transportation, Use, and Storage                  | 6         |
| Intended Conditions of Use                                                    | 7         |
| <b>SPYGLASS DS DIGITAL CONTROLLER COMPATIBILITY</b>                           | <b>7</b>  |
| <b>SETUP AND OPERATION</b>                                                    | <b>7</b>  |
| Front Panel Features                                                          | 7         |
| Rear Panel Features                                                           | 9         |
| Patient Environment                                                           | 10        |
| Note about Performance                                                        | 10        |
| Isolating the Controller from the Main Power Supply                           | 10        |
| Software                                                                      | 10        |
| Setup of the Controller                                                       | 10        |
| Starting the Controller                                                       | 12        |
| Powering On the Controller                                                    | 12        |
| Configure the Monitor for the Video Mode                                      | 13        |
| Adjusting the Brightness of the Video Image and Operating Light ON/OFF button | 14        |
| Recovering from a Controller Failure                                          | 15        |
| Considerations When Using a Laser or EHL Device                               | 15        |
| <b>COMPLETING A PROCEDURE</b>                                                 | <b>15</b> |
| Using the Controller In an Endoscopy Procedure                                | 15        |
| Shutting Down the Controller                                                  | 16        |
| Routine Inspection and Maintenance                                            | 16        |
| Cleaning and Disinfection                                                     | 16        |
| Disposal of the Product, Accessories, and Packing Materials                   | 16        |
| Post-Procedure                                                                | 16        |
| <b>INFORMATION TO BRIEF THE PATIENT</b>                                       | <b>16</b> |
| <b>TROUBLESHOOTING AND RESPONDING TO ERROR CODES</b>                          | <b>17</b> |
| Troubleshooting Chart                                                         | 17        |
| <b>WARRANTY</b>                                                               | <b>18</b> |
| <b>APPENDICES</b>                                                             | <b>18</b> |

## **Rx ONLY**

**Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician.

## **INTRODUCTION**

This user manual describes how to appropriately and safely use, maintain, and troubleshoot the SpyGlass™ DS Digital Controller (hereafter called the Controller).

The controller is used with a SpyScope™ DS Access and Delivery Catheter (M00546600) or SpyScope DS II Access and Delivery Catheter (M00546610). The SpyScope DS Access and Delivery Catheter and SpyScope DS II Access and Delivery Catheter are also called the SpyScope DS Catheter in this user manual.

## **Where to Get Help**

For technical support, ordering, service, and return authorization, contact Boston Scientific at 800-949-6708.

## **DEVICE DESCRIPTION**

The controller is an electronic device that:

- Receives video signals from a SpyScope DS Catheter,
- Processes the video signals, and
- Outputs video images to a video monitor.

The controller also generates and controls the light transmitted to the tip of the SpyScope DS Catheter to illuminate the area of interest within the anatomy. Buttons on the controller's front panel enable you to control the brightness level of the light.

To use the controller, connect it to a video monitor with a video cable and then connect a SpyScope DS Catheter to the controller. The controller provides direct visualization of the pancreatobiliary duct anatomy and enables exploratory and endotherapy procedures.

## **Contents**

- One (1) SpyGlass DS Digital Controller
- One (1) 3.1-meter North America power cord (M00546650) or one (1) 3.1-meter Brazil power cord (M0054665B0)
- One (1) 2.0-meter DVI cable
- One (1) 1.8-meter VGA cable

Ensure the package contains the components listed above.

## **Model Numbers**

|                                                                |            |
|----------------------------------------------------------------|------------|
| SpyGlass DS Digital Controller (with North America power cord) | M00546650  |
| SpyGlass DS Digital Controller (with Brazil power cord)        | M0054665B0 |

## **Specifications**

| <b>Electrical Specifications</b> |                              |
|----------------------------------|------------------------------|
| Input voltage                    | 100 to 240 VAC, 50/60 Hz     |
| Rated current                    | 1A-0.5A                      |
| Fuse rating                      | 250 V, 2A, Type F (F2AH250V) |

| Power Cord Specifications         |                                                                                                                                   |
|-----------------------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| 115 VAC<br>Domestic (US)          | Length—3.1 meters (10 feet)<br>Voltage Rating—125 VAC<br>Current Rating—10 amps<br>Connector Type—IEC 60320 C13                   |
| 230 VAC<br>International (Brazil) | Length—3.1 meters (10 feet)<br>Voltage Rating—250 VAC<br>Current Rating—10 amps<br>Connector Type—IEC 60320 C13<br>UL Marking—Yes |
| Physical Specifications           |                                                                                                                                   |
| Height                            | 11.5 cm (4.5 in)                                                                                                                  |
| Width                             | 33.0 cm (13.0 in)                                                                                                                 |
| Depth                             | 39.5 cm (15.5 in)                                                                                                                 |
| Weight (unpackaged)               | 6.8 kg (15 lb)                                                                                                                    |

### User Information

The controller and these instructions are intended for use by physicians trained in endoscopic pancreaticobiliary procedures, including endoscopic retrograde cholangio-pancreatography (ERCP).

A thorough understanding of the techniques, principles, clinical applications and risks associated with ERCP and cholangio-pancreatography procedures is suggested before using the controller as part of the SpyGlass DS and DS II Direct Visualization System (hereafter referred to as the SpyGlass DS System).

### INTENDED USE/INDICATIONS FOR USE

The SpyGlass DS and DS II Direct Visualization System is indicated for use in diagnostic and therapeutic applications during endoscopic procedures in the pancreatico-biliary system including the hepatic ducts. The SpyGlass DS and DS II Direct Visualization System comprises two components: the SpyScope DS Access and Delivery Catheter or SpyScope DS II Access and Delivery Catheter, and the SpyGlass DS Digital Controller.

The SpyScope DS Access and Delivery Catheter and SpyScope DS II Access and Delivery Catheter are intended to provide direct visualization and to guide both optical and accessory devices for diagnostic and therapeutic applications during endoscopic procedures in the pancreatico-biliary system including the hepatic ducts.

The SpyGlass DS Digital Controller is intended to provide illumination and receive, process, and output images from the SpyScope DS Access and Delivery Catheter or SpyScope DS II Access and Delivery Catheter for diagnostic and therapeutic applications during endoscopic procedures in the pancreaticobiliary system including the hepatic ducts.

### Clinical Benefit Statement

The SpyGlass DS Digital Controller when used with the SpyScope DS Access and Delivery Catheter or SpyScope DS II Access and Delivery Catheter provides direct optical visualization to access a target site and the ability to guide accessory devices within the pancreatico-biliary system including the hepatic ducts, allowing for delivery of diagnostic and therapeutic applications.

### CONTRAINDICATIONS

Contraindications associated with the use of this device include:

- Patients for whom ERCP is medically contraindicated.
- Contraindications specific to endoscopic pancreatico-biliary duct exploration and cannulation.

---

## WARNINGS

- Read this user manual, the SpyScope DS Catheter directions for use, and the monitor user manual before using the controller. Failure to follow any instructions or failure to heed any warnings or precautions may result in harm or injury to patient.
- Do not use the controller in the presence of flammable fluids and gases such as alcohol or oxygen. Doing so can result in fire and burns to the operator and patient.
- The LEDs in the catheter cable receptacle remain hot for a period of time after use. To prevent burns, do not insert fingers into catheter cable receptacle.
- Do not perform diagnostic or therapeutic procedures without a clear and adequate video display. Doing so can result in adverse events.
- Placing the controller where other electrical medical devices can degrade the video image can delay the procedure and result in adverse events. In addition, placing the controller where it can degrade the performance of other equipment in the endoscopy suite due to EMI emissions can delay the procedure or result in adverse events. To ensure the controller displays a clear and adequate video image and does not degrade the performance of other equipment, locate the controller as described in Table 6, Table 7, Table 8, and Table 9 in the Appendices. Verify operation in the endoscopy suite environment before starting a procedure. Follow ancillary equipment directions for use to locate ancillary equipment.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Using a controller without disinfecting the cabinet and front panel buttons can expose the operator to biohazardous materials. To prevent exposure to biohazardous materials, disinfect the chassis between use using the cleaning procedure described in "Cleaning and Disinfection."
- If the controller is connected to an improperly grounded power supply, electrical leakage can result in electrical shock to the user. To avoid risk of electrical shock, this equipment must only be connected to a supply mains with protective earth.
- Do not touch connecting devices for electrical connections between the different components (such as signal input and output connections for video signals, data exchange, control circuits, etc.) and the patient at the same time. Doing so can result in electric shock to the patient.
- If the controller experiences an unintended shutdown or lock up during a procedure, follow the procedure described in "Recovering from a Controller Failure." Failure to follow this recovery procedure after a controller failure can result in patient injury.
- No modification of this equipment is allowed.
- The use of accessories and cables other than those specified or supplied as spare parts from Boston Scientific may result in increased emissions or decreased immunity of the controller or SpyGlass DS System.
- Components added to the SpyGlass DS System by the user must be certified to the respective IEC standards (IEC 60601-1 for medical equipment, IEC 60950 for data processing equipment, and IEC 60065 for A/V equipment). In addition, the user must ensure the new configuration complies with the IEC 60601-1 system standard. Power cords must be compliant with the system.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the SpyGlass DS System, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

- The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

## PRECAUTIONS

- If you block the controller's ventilation outlet, the controller can overheat, resulting in a thermal shutdown or equipment damage. Leave at least 12.7 mm (0.5 in) between the controller back panel and other objects and 12.7 mm (0.5 in) of space between the side panels and other objects. Use the controller on a dedicated equipment cart to ensure proper ventilation.
- Spilling liquids on the controller can damage it or cause it to shut down. Do not place liquids above or near the controller.
- Opening the cabinet for repair purposes can damage the controller. The controller does not use operator-serviceable components. To prevent damage, do not access the controller cabinet.
- Connecting catheters other than a SpyScope DS Catheter to the controller can damage the controller. Only connect a SpyScope DS Catheter to the controller. See the "SpyGlass DS Digital Controller Compatibility" section.
- Locate the controller appropriately to avoid accidentally pulling cable connections, which can result in disconnection and loss of visualization.
- Before starting a procedure, ensure components such as the monitor and irrigation pump that support the SpyGlass DS System are present and operational. Starting a procedure without the supporting components present and operational can prolong the procedure.
- Do not use cleaners or disinfectant solutions that contain long-life surfactants. Doing so can leave conductive residues on the contacts of the catheter connector receptacle. The conductive residues can lead to malfunctions of the controller.
- Do not clean the LEDs located inside the catheter cable receptacle.
- Use of a cardiac defibrillator while a SpyScope DS Catheter remains in a patient can damage the controller. To prevent damage to the controller when using a defibrillator, remove the SpyScope DS Catheter before using the defibrillator.
- Do not insert a wet connector into the controller receptacle as poor video performance or damage to the controller may result.
- Applied parts of other electrical medical equipment in application with this equipment have to be type BF. Hence, only use the controller with the SpyScope DS Catheter.

## ADVERSE EVENTS

See directions for use for SpyScope DS Catheter for adverse event information.

## CONFORMANCE TO STANDARDS

### Essential Performance Statement

The controller does not have essential performance as defined by IEC 60601-1 and IEC 60601-2-18.

## HOW SUPPLIED

Device supplied non-sterile. Inspect the components for damage. Do not use a component if it appears damaged. Do not use the components if the package labeling is incomplete or illegible.

## Handling and Storage

### Environmental Limits During Transportation, Use, and Storage

The controller can be stored between procedures on the endoscopy tower or on a dedicated equipment cart (Table 1).

**Table 1. Transportation, use, and storage environmental limits**

| Environmental Limits during Use                   |                           |
|---------------------------------------------------|---------------------------|
| Ambient temperature (°C)                          | 10 to 35                  |
| Relative humidity (%)                             | 30 to 85 (non-condensing) |
| Atmospheric pressure (hPa)                        | 700 to 1060               |
| Environmental Limits during Transport and Storage |                           |
| Ambient temperature (°C)                          | -40 to 70                 |
| Relative humidity (%)                             | 10 to 90 (non-condensing) |
| Atmospheric pressure (hPa)                        | 500 to 1060               |

### Intended Conditions of Use

The controller should be used in a hospital.

### SPYGLASS DS DIGITAL CONTROLLER COMPATIBILITY

All attached equipment must comply with its applicable IEC standards.

The controller is compatible with:

- SpyScope DS Access and Delivery Catheter (M00546600).
- SpyScope DS II Access and Delivery Catheter (M00546610) - compatible with Controllers with software revision 2.2 only.

The controller incorporates standard video signal outputs (i.e. S-Video, DVI, SXGA) to allow connection to a monitor. Boston Scientific has tested the following monitors to ensure compatibility with the controller and IEC standards:


- Sony LMD1951MD monitor
- Sony LMD1950MD monitor

An irrigation pump is required to perform procedures with the SpyGlass DS System. Boston Scientific has tested the SpyGlass Irrigation Pump to ensure it is compatible with the SpyGlass DS System.

### SETUP AND OPERATION


#### Front Panel Features

Figure 1, Table 2, and Table 3 illustrate and describe the controller's front panel features.

**Figure 1. Illustration of front panel features**

**Table 2. Description of front panel features**

|   | Feature                                                                             | Description                                                                                                                                                                                                                                                                                                                                                                                                       |
|---|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 |    | <b>Power Button:</b> When the controller is connected to a power source, pushing the power button alternates between powering up and powering down the controller. The power button illuminates green when the controller is powered up.                                                                                                                                                                          |
| 2 |    | <b>Connection Cable Receptacle:</b> You plug the connection cable from the SpyScope DS Catheter into the receptacle during setup of the system.                                                                                                                                                                                                                                                                   |
| 3 |    | Indicates you must use BF applied parts.                                                                                                                                                                                                                                                                                                                                                                          |
| 4 |    | <b>Light On/Off Button:</b> Pressing this button alternately turns the light off and on. Light On/Off Button—When the controller is on, each press of the Light On/Off button toggles the light between on and off. When the light is on, the button illuminates blue. When the light is off, the button illuminates white.                                                                                       |
| 5 |   | <b>Light Brightness Control Buttons:</b> When the light is on, the light brightness control buttons are illuminated blue. Pressing the  button increases the intensity of the light. Pressing the  button decreases the intensity of the light. |
| 6 |  | <b>Light Brightness Indicator:</b> The bars of this indicator illuminate to indicate the brightness of the light. There are five levels of brightness (Table 3).                                                                                                                                                                                                                                                  |

**Table 3. Interpreting the light brightness indicator**

| Illumination State of Light Brightness Indicator                                    | Setting Description |
|-------------------------------------------------------------------------------------|---------------------|
|  | Low                 |
|  | Med-Low             |
|  | Med                 |
|  | Med-High            |
|  | High                |

## Rear Panel Features

Figure 2 and Table 4 illustrate and describe the controller's rear panel features.


Figure 2. Illustration of back panel features

Table 4. Description of back panel features


|   | Feature                                                                             | Description                                                                                                                                                                                                                                                                           |
|---|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 |   | USB port access cover (tool required).                                                                                                                                                                                                                                                |
| 2 |  | <b>DVI Output Connector:</b> Video output connector for compatible DVI-capable monitors.                                                                                                                                                                                              |
| 3 |  | Symbol indicating video outputs.                                                                                                                                                                                                                                                      |
| 4 |  | <b>VGA Output Connector:</b> Video output connector for compatible SXGA capable monitors.                                                                                                                                                                                             |
| 5 |  | <b>Two S-video Output Connectors:</b> Video output connectors for compatible S-video-capable monitors or accessories. [A user-supplied S-Video cable is required to connect the SpyGlass DS Controller to a monitor. The recommended S-Video cable is TecNec part number SV4-SV4-15.] |
| 6 |  | <b>PAL/NTSC Switch:</b> The switch enables operation of the controller with PAL or NTSC video formats when using an S-video cable. PAL is the Phase Alternating Line format. NTSC is the National Television System Committee format.                                                 |
| 7 |  | <b>Cabinet Ventilation Outlet:</b> The cabinet ventilation outlet enables the cooling fan to exhaust hot air from the cabinet to maintain the proper operating temperature inside the cabinet.                                                                                        |

|    | Feature                                                                           | Description                                                                                                                                                                                                                                          |
|----|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 8  |  | Symbol for potential equalization conductor.                                                                                                                                                                                                         |
| 9  |  | <b>Potential Equalization Conductor:</b> Provides a means of securely linking the potential equalization conductor of the controller to ground.                                                                                                      |
| 10 |  | <b>Fuse Holder Module and Power Cord Connector:</b> The fuse holder module provides access to the controller electrical fuses. The power cord connector accepts the AC power cord that connects to an AC Main or the SpyGlass Isolation Transformer. |
| 11 |  | <b>Label:</b> Provides regulatory and manufacturing information.                                                                                                                                                                                     |

### Patient Environment

The controller is patient equipment and can be used in the patient environment (Figure 3).

Place the controller to allow easy access to the power cord in the event you must quickly disconnect the controller from the main power source

**Figure 3. Typical positioning of equipment, patient, and operator**

### Note about Performance

To ensure the controller performs within its design specifications, position it following the guidelines presented in the Appendices.

### Isolating the Controller from the Main Power Supply

To isolate the controller from the main power supply, disconnect the power cord from the power supply receptacle or from the power cord receptacle on the controller.

### Software

The revision level of the software installed on your controller is displayed on the cable connect screen.

### Setup of the Controller

All attached equipment must comply with its applicable IEC standards.

Upon receiving the controller, complete these initial setup steps:

1. Inspect the controller and its components to ensure they are not damaged.
2. Clean the controller following the instructions in the "Cleaning and Disinfection" section.
3. Place the controller and the video monitor on an equipment cart such as the SpyGlass Components Travel Cart or on an existing endoscopy tower in the procedure room. See the "SpyGlass DS Digital Controller Compatibility" section for additional information if needed.
4. Set the PAL/NTSC video standard switch on the back of the controller as follows:
  - a. Use of DVI or VGA cable—No PAL/NTSC switch adjustment is required.
  - b. Use of S-video cable—Set the PAL/NTSC video mode switch on the back of the controller to match the video mode of the monitor (Figure 4).
5. Use the DVI or VGA cable (or a user-supplied S-Video cable) to connect one of the video input connectors on the video monitor to one of the video output connectors on the controller.

---

**Note:** Connecting the monitor to the controller using the DVI cable provides the highest quality video image.

---

6. Connect the controller and video monitor power cords to electrical receptacles or an optional SpyGlass Isolation Transformer.
7. Power on and set the monitor video input based on the cable selected in Step 5. Follow the steps below to set the monitor's video input.

**Sony LMD 1951 Monitor:**

---

**Note:** All references to buttons refer to the buttons along the front panel of the monitor.

---

1. Press **Control** to unlock the front panel buttons
  - a. Select the **DVI** button for DVI operation
  - b. Select the **HD15** for SXGA operation
  - c. Select the **Y/C** button for S-Video operation

**Sony LMD 1950 Monitor:**

---

**Note:** All references to buttons refer to the buttons along the front panel of the monitor.

---

- a. Push the **Input** up arrow to select **HD15** for SXGA operation
  - b. Push the **Input** up arrow to select **DVI** for DVI operation
  - c. Push the **Input** up arrow to select **Y/C** for S-Video operation
8. Power on the controller using the steps in the "Starting the Controller" section.

---

**Note:** Controller must be powered on in order to configure monitor, however SpyScope DS Catheter does not need to be connected.

---


9. Configure the monitor for use in DVI, SXGA, or S-video mode following the steps in the "Configure the Monitor for the Video Mode" section.

---

**Note:** Monitor must be configured before beginning a procedure.

---

**Figure 4. Back panel connection points**

### Starting the Controller

Follow the steps below to start the controller. You can start the controller with or without the catheter connection cable plugged into the controller.

### Powering On the Controller

1. Press the power button to start the controller. The power button will illuminate green, and the controller will begin a self-test and boot sequence. The monitor displays the **boot screen** (a) followed by the **transition screen** (b). If monitor does not display the boot screen, consult the trouble shooting section.

2. Upon a successful boot up, the monitor displays the **cable connect** screen (c) showing the software version (d). Connect the catheter connection cable to the controller until it locks into place. Pull back on the cable connector to ensure the connector is fully seated in the controller.


If the SpyScope DS Catheter is already connected, this screen will not appear, proceed to step 3.

---


**Note:** If you do not want the Boston Scientific watermark (e) to appear on the screen, simultaneously press the  and  buttons for at least 3 seconds while the cable connect screen is displayed on the monitor. To redisplay the watermark, press the  and  buttons for at least 3 seconds.

---


3. The controller briefly displays the **catheter loading** screen (f).

4. Shortly after the catheter loading screen appears, a live video image will appear on the live video screen (g). If the live image does not appear, consult the trouble shooting section.

### Configure the Monitor for the Video Mode

Configure the monitor according to the monitor's user guide. Recommendations for the Sony LMD1951MD and LMD1950MD monitors are as follows:

**Option 1—(Recommended)**—Configure monitor for SXGA or DVI mode.

**Option 2**—Configure monitor for S-Video mode.

### Option 1—Configure monitor for SXGA or DVI mode.

---

**Note:** All references to buttons refer to the buttons along the front panel of the monitor.

---

## Configure the LMD1951MD Monitor for SXGA or DVI

Recall the default monitor setting by:

- a. Pressing **User Mem** to open the **User Memory Screen**.
- b. Highlight **Default** in the **User Memory Screen** window.
- c. Press **Enter**.

## Configure the LMD1950MD Monitor for SXGA or DVI

Recall the default monitor setting by:

- a. Pressing **User Memory** to open the **User Memory Screen**.
- b. Highlight **Default** in the **User Memory Screen** window.
- c. Press **Enter**.

## Option 2—Configure monitor for S-Video mode.

---

**Note:** All references to buttons refer to the buttons along the front panel of the monitor.

---

## Configure the LMD1951MD Monitor for S-Video

1. Recall the default monitor setting by:
  - a. Pressing **User Mem** to open the **User Memory Screen**.
  - b. Highlight **Default** in the **User Memory Screen** window.
  - c. Press **Enter**.
2. Press **Menu** to open the status window.
3. Press the menu '–' button multiple times until the **User Config** window is displayed.
4. Press **Enter** to highlight the **System Setting** line in yellow.
5. Press **Enter** to open the **User Config** window.
6. Press the menu '–' button multiple times to highlight the **I/P Mode** line in yellow.
7. Press **Enter** to highlight the **I/P Mode**.
8. Press the menu '–' button until the **I/P Mode** is set to **Field Merge**.
9. Press **Enter**.
10. Press **Menu** to close the screen menus.

## Configure the LMD1950MD Monitor for S-Video


Recall the default monitor setting by:

- a. Pressing **User Memory** to open the **User Memory Screen**.
- b. Highlight **Default** in the **User Memory Screen** window.
- c. Press **Enter**.

## Adjusting the Brightness of the Video Image and Operating Light ON/OFF button

To increase the brightness of the video image as displayed on the **live video** screen, press the  button. To decrease the brightness of the video image, press the  button. The image brightness scale provides visual feedback of the relative brightness of the light source.

Pressing the light On/Off button toggles the controller light source between ON and OFF state. To operate the light On/Off button the controller must be powered on and a SpyScope DS Catheter must be connected to the controller. Button illuminates white when the light is off and blue when it is on.

## Recovering from a Controller Failure

There are two possible failure modes: (1) unintended controller shutdown, and (2) controller lock up.

### Recovery from Unintended Shutdown

Follow these steps to recover from an unintended controller shutdown:

1. Remove any accessory from the SpyScope DS Catheter.
2. Remove the SpyScope DS Catheter from the patient.
3. Unplug the SpyScope DS Catheter from the controller.
4. Press the controller's **Power** button to restart the controller.
5. If the controller fails to start, contact Boston Scientific.

### Recovery from a Controller Lock Up

Follow these steps to recover from a controller lock up:

1. Remove any accessory from the SpyScope DS Catheter.
2. Remove the SpyScope DS Catheter from the patient.
3. Unplug the SpyScope DS Catheter from the controller.
4. Press and hold the power button until the controller shuts down.
5. Press the controller's **Power** button to restart the controller.
6. If the controller fails to start, contact Boston Scientific.

## Considerations When Using a Laser or EHL Device

When activating energy from a laser or electrohydraulic lithotripsy (EHL) generator during lithotripsy, you can expect to see a bright flash on the video monitor and possibly a momentary disruption in the video quality. This reaction is normal and is not indicative of a defect or operating problem with the controller or SpyScope DS Catheter.

## COMPLETING A PROCEDURE

The procedure below describes the procedural use of the controller and assumes you have received, inspected, assembled, and tested the controller following instructions in the "Setup and Operation" section.

### Using the Controller In an Endoscopy Procedure

Using the controller includes the following steps:

1. Clean the controller as directed in the "Cleaning and Disinfection" section.
2. Power up the controller.
3. Configure the monitor.
4. Connect a SpyScope DS Catheter to the front panel connector.
5. Complete the procedure as outlined in SpyScope DS Catheter directions for use.

## Shutting Down the Controller

To shut down the controller at the end of a procedure or during a procedure, follow these steps:

1. Remove any accessory and the SpyScope DS Catheter from the patient following the instructions in the SpyScope DS Catheter directions for use.
2. Disconnect the catheter cable from the front of the controller by pushing down on the cable connector locking tab and pulling it out of the receptacle.
3. Power down the controller by pressing the Power button. The power button indicator light shuts off indicating power has been cut to the controller.
4. If you are shutting the controller down after finishing a procedure, dispose of the SpyScope DS Catheter as described in the SpyScope DS Catheter directions for use. Then clean the controller as directed in the "Cleaning and Disinfection" section.

## Routine Inspection and Maintenance

The controller does not require routine maintenance and calibration because the self-test mode, activated automatically upon turning on the power to the controller, verifies correct functioning of the controller. Periodically check the power cord assembly for damage to the insulation or connectors. In the event that the controller requires repair or replacement, contact Boston Scientific.

## Cleaning and Disinfection

Disconnect the power cord before cleaning or disinfecting the unit. Use a 15 to 70 percent isopropyl alcohol in purified water solution and a cloth to clean the controller enclosure, front panel, and power cable. Do not allow fluids to enter the enclosure, power cable connections, catheter cable receptacle, or component/accessory connections. Do not attempt to clean the unit while it is plugged into an electrical outlet.

## Disposal of the Product, Accessories, and Packing Materials

1. Clean and disinfect all external and accessible surfaces using the cleaning and disinfection instructions included in the user manual. Include any common detachable cables (power cord and video cables.)
2. Use an electronics recycling facility to dispose of the Controller. When placing the unit into an electronics recycling stream, notify the receiver of the presence of lithium ion coin batteries.
  - a. Use of recycling service suppliers familiar with medical electrical equipment is recommended, but not required.

Do not dispose of by incineration, burial or placement into a common waste stream.

Note: for users in California, USA. Perchlorate Material—special handling may apply. Battery in controller contains Perchlorate. See [www.dtsc.ca.gov/hazardouswaste/perchlorate](http://www.dtsc.ca.gov/hazardouswaste/perchlorate).

For further information, contact Boston Scientific.

## Post-Procedure

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

## INFORMATION TO BRIEF THE PATIENT

Fully explain to the patient the potential benefits and risks of the endoscopy and endoscopic treatment as well as any examination/treatment methods that can be performed in its place, and perform the endoscopy and endoscopic treatment only after obtaining the consent of the patient.

Even after starting the endoscopy and endoscopic treatment, continue to evaluate the potential benefits and risks, and immediately stop the endoscopy/treatment and take proper measures if the risks to the patient become greater than the potential benefits.

## TROUBLESHOOTING AND RESPONDING TO ERROR CODES

### Troubleshooting Chart

Most operating problems are easily resolved. If the controller does not operate as expected, try to resolve the problem with this troubleshooting chart before contacting Boston Scientific for technical support (Table 5).

**Table 5. Controller troubleshooting chart**

| Symptom                                                                                                          | Possible Cause                                             | Corrective Action                                                                                                                                                      |
|------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Controller does not start or power button is not illuminated green when pressed                                  | Power cord not connected or not connected tightly          | Ensure both ends of the power cord are firmly connected to their connection points.                                                                                    |
|                                                                                                                  | No electrical power at the receptacle                      | Check the receptacle's circuit breaker to ensure it is not tripped.                                                                                                    |
|                                                                                                                  |                                                            | Confirm the receptacle has power by operating another electrical device from the receptacle.                                                                           |
|                                                                                                                  | Boot problem                                               | Reboot the controller by powering it down and restarting it. If the problem reoccurs, call Boston Scientific for support.                                              |
|                                                                                                                  | Defective power cord                                       | Replace power cord.                                                                                                                                                    |
|                                                                                                                  | Controller fuse blown                                      | Call Boston Scientific for support.                                                                                                                                    |
|                                                                                                                  | Controller damaged                                         | Call Boston Scientific for support.                                                                                                                                    |
| Catheter connection cable is plugged into the controller, but no video image is displayed.                       | Video monitor is not powered up or not connected           | Power up the video monitor. Check that the video cable is properly connected to the monitor and controller. Verify that the monitor is set to the correct video input. |
|                                                                                                                  | Catheter cable is not properly connected to controller     | Ensure the plug on the catheter connection cable is firmly inserted into the connector with the locking tab facing up. Ensure connection terminals are dry and clean.  |
|                                                                                                                  | SpyScope DS Catheter is broken or defective                | Replace SpyScope DS Catheter.                                                                                                                                          |
|                                                                                                                  | Illumination light is off                                  | Press illumination light button to turn on light. (The light On/Off button should be illuminated blue.)                                                                |
| Catheter connection cable is plugged into the controller, but <b>cable connection</b> screen is still displayed. | Catheter cable is not properly connected to controller     | Ensure the plug on the catheter connection cable is firmly inserted into the connector with the locking tab facing up. Ensure connection terminals are dry and clean.  |
|                                                                                                                  | SpyScope DS Catheter is broken or defective                | Replace SpyScope DS Catheter.                                                                                                                                          |
| Video image is too dark.                                                                                         | Brightness setting is too low                              | Adjust brightness using brightness control buttons.                                                                                                                    |
|                                                                                                                  | Debris covering the distal tip of the SpyScope DS Catheter | Clean the distal tip by irrigating, or removing & wiping it with 15 to 70 percent isopropyl alcohol in purified water solution using a cotton swab.                    |
|                                                                                                                  | SpyScope DS Catheter fiber optic is damaged                | Replace SpyScope DS Catheter.                                                                                                                                          |
| Video image is too bright                                                                                        | Brightness setting is too high                             | Adjust brightness using brightness control buttons.                                                                                                                    |
|                                                                                                                  | Video cable connect to monitor output                      | Ensure video cable is connected into monitor video input and not output.                                                                                               |

| Symptom                                                                                     | Possible Cause                                                                  | Corrective Action                                                                                                                                                                                                                                            |
|---------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Video image is blurry, scrambled, inadequate, or distorted                                  | Controller is too close to other electrical medical equipment                   | Ensure controller is placed as directed in Table 9 of Appendix 4.<br>Shutdown other medical electrical equipment to determine which equipment is causing the symptom. Properly place other electrical medical devices according to their directions for use. |
|                                                                                             | Video cable not completely connected to either monitor or controller            | Ensure video cable is completely connected to video monitor and controller.                                                                                                                                                                                  |
| Video image is scrambled or displays unnatural colors                                       | PAL/NTSC switch is set improperly                                               | Determine which video mode standard your monitor uses and set the PAL/NTSC switch accordingly using the recessed switch on the rear panel. Once switch is set correctly power on the controller.                                                             |
| Video image is blurry, distorted, or otherwise not acceptable                               | Debris covering the distal tip of the SpyScope DS Catheter                      | Clean the distal tip by irrigating or removing & wiping it with 15 to 70 percent isopropyl alcohol in purified water solution using a cotton swab.                                                                                                           |
|                                                                                             | Video monitor is not compatible with controller                                 | Replace monitor with compatible monitor.                                                                                                                                                                                                                     |
| Chassis is warm-to-hot to the touch                                                         | Chassis ventilation outlet is blocked with debris or too close to other objects | Reposition controller to provide more clearance for ventilation. If not resolved, contact Boston Scientific for support.                                                                                                                                     |
| If any of these problems persist, please contact BSC for repair or replacement information. |                                                                                 |                                                                                                                                                                                                                                                              |

## WARRANTY

For device warranty information, visit ([www.bostonscientific.com/warranty](http://www.bostonscientific.com/warranty)).

EU Importer: Boston Scientific International B.V., Vestastraat 6, 6468 EX Kerkrade, The Netherlands

SpyGlass and SpyScope are trademarks of Boston Scientific Corporation or its affiliates.

All other trademarks are the property of their respective owners.

## APPENDICES

The appendices include:

- Appendix 1—Guidance and Manufacturer's Declaration-Electromagnetic Emissions
- Appendix 2—Guidance and Manufacturer's Declaration-Electromagnetic Immunity
- Appendix 3—Guidance and Manufacturer's Declaration— electromagnetic IMMUNITY – for ME Equipment and ME Systems that are not LIFE-SUPPORTING
- Appendix 4—Recommended Separation Distance Between Portable and Mobile RF Communications Equipment and the SpyGlass DS Digital Controller
- Appendix 5—Medical Design Criteria and Specifications

## Appendix 1—Guidance and Manufacturer's Declaration-Electromagnetic Emissions

The SpyGlass DS Digital Controller is intended for use in the electromagnetic environment specified below (Table 6). The customer or the user of the SpyGlass DS Digital Controller should assure that it is used in such an environment.

**Table 6. Guidance and manufacturer's declaration - electromagnetic emissions**

| Emissions Test                                          | Compliance     | Electromagnetic Environment—Guidance                                                                                                                                                                                                                                                                                                                                                                                                            |
|---------------------------------------------------------|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| RF Emissions<br>CISPR 11                                | Group 1        | The SpyGlass DS Digital Controller uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.<br><br>The SpyGlass DS Digital Controller is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| RF Emissions<br>CISPR 11                                | Class A        |                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Harmonic Emissions<br>IEC 61000-3-2                     | Not Applicable |                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Voltage Fluctuations/Flicker<br>Emissions IEC 61000-3-3 | Not Applicable |                                                                                                                                                                                                                                                                                                                                                                                                                                                 |

## Appendix 2—Guidance and Manufacturer's Declaration-Electromagnetic Immunity

The SpyGlass DS Digital Controller is intended for use in the electromagnetic environment specified below (Table 7). The customer or the user of the SpyGlass DS Digital Controller should assure that it is used in such an environment.

**Table 7. Guidance and manufacturer's declaration—electromagnetic immunity IEC 60601-1-2 (2nd edition and 3rd edition)**

| Immunity Test                                                                                          | IEC 60601 Test Level                                                                                                                                                                | Compliance Level                                                                                                                                                                     | Electromagnetic Environment (Guidance)                                                                                                                                                                                                                                                                                       |
|--------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Electrostatic discharge (ESD)<br>IEC 61000-4-2                                                         | ± 6 kV contact<br>± 8 kV/V air                                                                                                                                                      | ± 6 kV contact<br>± 8 kV/V air                                                                                                                                                       | Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.                                                                                                                                                                               |
| Electrical fast transient / burst<br>IEC 61000-4-4                                                     | ± 2 kV for power supply lines<br>± 1 kV for input/output lines                                                                                                                      | ± 2 kV for power supply lines<br>± 1 kV for input/output lines                                                                                                                       | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                                                                          |
| Surge<br>IEC 61000-4-5                                                                                 | ± 1 kV differential mode<br>± 2 kV common mode                                                                                                                                      | ± 1 kV differential mode<br>± 2 kV common mode                                                                                                                                       | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                                                                          |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br>IEC 61000-4-11 | < 5% UT<br>(> 95% dip in UT for 0.5 cycle<br><br>40% U<br>(60% dip in UT) for 5 cycles<br><br>70% UT<br>(30% dip in UT) for 25 cycles<br><br>< 5% UT<br>(> 95% dip in UT) for 5 sec | < 5% UT<br>(> 95% dip in UT) for 0.5 cycle<br><br>40% U<br>(60% dip in UT) for 5 cycles<br><br>70% UT<br>(30% dip in UT) for 25 cycles<br><br>< 5% UT<br>(> 95% dip in UT) for 5 sec | Mains power quality should be that of a typical commercial or hospital environment. If the user of the SpyGlass DS Digital Controller requires continued operation during power mains interruptions, it is recommended that the SpyGlass DS Digital Controller be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60 Hz) magnetic field<br>IEC 61000-4-8                                             | 3 A/m                                                                                                                                                                               | 3 A/m                                                                                                                                                                                | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.                                                                                                                                                                                    |
| Note — UT is the a.c. mains voltage prior to application of the test level.                            |                                                                                                                                                                                     |                                                                                                                                                                                      |                                                                                                                                                                                                                                                                                                                              |

| Immunity Test                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          | IEC 60601 Test Level            | Compliance Level | Electromagnetic Environment (Guidance)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------|------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Conducted RF<br>IEC 61000-4-6                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          | 3 Vrms<br><br>150 kHz to 80 MHz | 3 Vrms           | Portable and mobile RF communications equipment should be used no closer to any part of the SpyGlass DS Digital Controller, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.<br><br><b>Recommended separation distance</b><br><br>$d=1.2\sqrt{P}$<br><br>$d=1.2\sqrt{P}$ 80 MHz to 800 MHz<br><br>$d=2.3\sqrt{P}$ 800 MHz to 2.5 GHz<br><br>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup> , should be less than the compliance level in each frequency range <sup>b</sup> . Interference may occur in the vicinity of equipment marked with the following Symbol:<br><br> |
| Radiated RF<br>IEC 61000-4-3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           | 3 V/m<br><br>80 MHz to 2.5 GHz  | 3 V/m            |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Note 1 — At 80 MHz and 800 MHz, the higher frequency range applies.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                 |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Note 2 — These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |                                 |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the SpyGlass DS Digital Controller is used exceeds the applicable RF compliance level above, the SpyGlass DS Digital Controller should be observed to verify normal operation.<br>If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the SpyGlass DS Digital Controller.<br>b) Over the frequency range 150 kHz to 80 MHz, the field strengths should be less than 3 V/m. |                                 |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |

### Appendix 3—Guidance and Manufacturer's Declaration-Electromagnetic Immunity

The SpyGlass DS Digital Controller is intended for use in the electromagnetic environment specified below (Table 8). The customer or the user of the SpyGlass DS Digital Controller should assure that it is used in such an environment.

**Table 8. Guidance and manufacturer's declaration—electromagnetic immunity IEC 60601-1-2 (4th edition)**

| Immunity Test                                                                  | IEC 60601<br>Test Level/Compliance Level                                  | Electromagnetic Environment<br>(Guidance)                                                                                                                                                                                                                                                                                    |
|--------------------------------------------------------------------------------|---------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Electrostatic discharge (ESD)<br>IEC 61000-4-2                                 | ± 8 kV contact<br>± 15 kV air                                             | Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.                                                                                                                                                                               |
| Electrical fast transient/burst<br>IEC 61000-4-4                               | ± 2 kV for power supply lines<br>± 1 kV for input/output lines            | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                                                                          |
| Surge<br>IEC 61000-4-5                                                         | ± 1 kV differential mode<br>± 2 kV common mode                            | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                                                                          |
| Voltage dips<br>IEC 61000-4-11                                                 | 0% $U_T$ ; 0.5 cycle<br>At 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° | Mains power quality should be that of a typical commercial or hospital environment. If the user of the SpyGlass DS Digital Controller requires continued operation during power mains interruptions, it is recommended that the SpyGlass DS Digital Controller be powered from an uninterruptible power supply or a battery. |
|                                                                                | 0% $U_T$ ; 1 cycle and<br>70% $U_T$ ; 25/30 cycles<br>Single phase: at 0° |                                                                                                                                                                                                                                                                                                                              |
| Voltage interruptions<br>IEC 61000-4-11                                        | 0% $U_T$ ; 250/300 cycle                                                  | Mains power quality should be that of a typical commercial or hospital environment. If the user of the SpyGlass DS Digital Controller requires continued operation during power mains interruptions, it is recommended that the SpyGlass DS Digital Controller be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60Hz) magnetic field<br>IEC 61000-4-8                      | 30 A/m                                                                    | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.                                                                                                                                                                                    |
| Note — $U_T$ is the a.c. mains voltage prior to application of the test level. |                                                                           |                                                                                                                                                                                                                                                                                                                              |

| Immunity Test                                                               | IEC 60601<br>Test Level/Compliance Level                                                   |             |                                       |                           | Electromagnetic Environment<br>(Guidance)                                                                                                                                                                                                                           |
|-----------------------------------------------------------------------------|--------------------------------------------------------------------------------------------|-------------|---------------------------------------|---------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Conducted RF<br>IEC 61000-4-6                                               | 3 V<br>0.15 MHz -80 MHz<br>6 V in ISM bands between<br>0.15 MHz -80 MHz<br>80% AM at 1 kHz |             |                                       |                           | Portable and mobile RF communications equipment should be used no closer to any part of the SpyGlass DS Digital Controller, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. |
| Radiated RF<br>IEC 61000-4-3                                                | 3 V/m<br>80 MHz to 2.7 GHz<br>80% AM at 1 kHz                                              |             |                                       |                           |                                                                                                                                                                                                                                                                     |
| Proximity fields from RF wireless communications equipment<br>IEC 61000-4-3 | Test Frequency (MHz)                                                                       | Band (MHz)  | Modulation                            | Immunity Test Level (V/m) | <b>Recommended separation distance</b><br><br>d=2√P 80 MHz to 2.7 GHz                                                                                                                                                                                               |
|                                                                             | 385                                                                                        | 380 - 390   | Pulse modulation<br>18 Hz             | 27                        |                                                                                                                                                                                                                                                                     |
|                                                                             | 450                                                                                        | 430 - 470   | FM<br>± 5 kHz deviation<br>1 kHz sine | 28                        | Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).                                                                                    |
|                                                                             | 710                                                                                        | 704 - 787   | Pulse modulation<br>217 Hz            | 9                         |                                                                                                                                                                                                                                                                     |
|                                                                             | 745                                                                                        |             |                                       |                           |                                                                                                                                                                                                                                                                     |
|                                                                             | 780                                                                                        |             |                                       |                           |                                                                                                                                                                                                                                                                     |
|                                                                             | 810                                                                                        | 800 - 960   | Pulse modulation<br>18 Hz             | 28                        | Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup> , should be less than the compliance level in each frequency range <sup>b</sup> .                                                                          |
|                                                                             | 870                                                                                        |             |                                       |                           |                                                                                                                                                                                                                                                                     |
|                                                                             | 930                                                                                        |             |                                       |                           |                                                                                                                                                                                                                                                                     |
|                                                                             | 1720                                                                                       | 1700 - 1990 | Pulse modulation<br>217 Hz            | 28                        | Interference may occur in the vicinity of equipment marked with the following Symbol:                                                                                                                                                                               |
|                                                                             | 1845                                                                                       |             |                                       |                           |                                                                                                                                                                                                                                                                     |
|                                                                             | 1970                                                                                       |             |                                       |                           |                                                                                                                                                                                                                                                                     |
|                                                                             | 2450                                                                                       | 2400 - 2570 | Pulse modulation<br>217 Hz            | 28                        |                                                                                                                                                                                |
|                                                                             | 5240                                                                                       | 5100 - 5800 | Pulse modulation<br>217 Hz            | 9                         |                                                                                                                                                                                                                                                                     |
| 5500                                                                        |                                                                                            |             |                                       |                           |                                                                                                                                                                                                                                                                     |
| 5785                                                                        |                                                                                            |             |                                       |                           |                                                                                                                                                                                                                                                                     |

Note — These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the SpyGlass DS Digital Controller is used exceeds the applicable RF compliance level above, the SpyGlass DS Digital Controller should be observed to verify normal operation.

If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the SpyGlass DS Digital Controller.

b) Over the frequency range 150 kHz to 80 MHz, the field strengths should be less than 3 V/m.

## Appendix 4—Recommended Separation Distance Between Portable and Mobile RF Communications Equipment and the SpyGlass DS Digital Controller

The SpyGlass DS Digital Controller is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the SpyGlass DS Digital Controller can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the SpyGlass DS Digital Controller as recommended below, according to the maximum output power of the communications equipment (Table 9).

**Table 9. Recommended separation distances between portable and mobile RF communications equipment and the SpyGlass DS Digital Controller IEC 60601-1-2 (2nd edition and 3rd edition)**

| Max Output Power (Watts)                                                                                                                                                                                                                                                                                        | Separation (m)<br>150 kHz to 80 MHz<br>$D=1.2(\text{Sqrt } P)$ | Separation (m)<br>80 MHz to 800 MHz<br>$D=1.2(\text{Sqrt } P)$ | Separation (m)<br>800 MHz to 2.5 GHz<br>$D=2.3(\text{Sqrt } P)$ |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------|----------------------------------------------------------------|-----------------------------------------------------------------|
| 0.01                                                                                                                                                                                                                                                                                                            | .12                                                            | .12                                                            | .23                                                             |
| 0.1                                                                                                                                                                                                                                                                                                             | .38                                                            | .38                                                            | .73                                                             |
| 1                                                                                                                                                                                                                                                                                                               | 1.2                                                            | 1.2                                                            | 2.3                                                             |
| 10                                                                                                                                                                                                                                                                                                              | 3.8                                                            | 3.8                                                            | 7.3                                                             |
| 100                                                                                                                                                                                                                                                                                                             | 12                                                             | 12                                                             | 23                                                              |
| For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power of the transmitter in watts (W) according to the manufacturer. |                                                                |                                                                |                                                                 |
| Note 1—At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.                                                                                                                                                                                                                   |                                                                |                                                                |                                                                 |
| Note 2—These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.                                                                                                                                              |                                                                |                                                                |                                                                 |

**Table 10. Recommended separation distances between portable and mobile RF communications equipment and the SpyGlass DS Digital Controller IEC 60601-1-2 (4th edition)**

| Max Output Power (Watts)                                                                                                                                                                                                                                                                                        | Separation Distance $D=2\sqrt{P}$ |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------|
| 0.01                                                                                                                                                                                                                                                                                                            | 0.2                               |
| 0.1                                                                                                                                                                                                                                                                                                             | 0.6                               |
| 1                                                                                                                                                                                                                                                                                                               | 2.0                               |
| 10                                                                                                                                                                                                                                                                                                              | 6.3                               |
| 100                                                                                                                                                                                                                                                                                                             | 20.0                              |
| For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power of the transmitter in watts (W) according to the manufacturer. |                                   |
| Note — These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.                                                                                                                                              |                                   |

## Appendix 5—Medical Design Criteria and Specifications

The SpyGlass DS Digital Controller meets medical design criteria and specifications (Table 11):

**Table 11. Medical design criteria and specification**

| <b>Design Criteria for SpyGlass DS Digital Controller</b> | <b>Specification</b>                        |
|-----------------------------------------------------------|---------------------------------------------|
| Type of protection against electric shock                 | Class 1 equipment                           |
| Degree of protection against electric shock               | Type BF applied part (SpyScope DS Catheter) |
| Degree of protection against ingress of materials         | IP20                                        |
| Mode of operation                                         | Continuous operation                        |
| Installation and use                                      | Portable equipment                          |

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All rights reserved.
© 2020 Boston Scientific Corporation or its affiliates.