

FibroScan[®]

compact 530 by echosens



Multi department use

- Fully featured mobility
- Non-invasive & quantitative liver exam

- Easy to handle
- VCTE™ liver fibrosis & cirrhosis assessment
- User-friendly tactile interface
- Battery operated device



- Dedicated roll stand
- Robust transport case
- 2 probes holders and connectors
- CAP™ Steatosis assessment

BENEFITS

- Non-invasive
 - Precise & reproducible
 - Standardized procedure
 - Easy to use
 - Can be performed by any trained operator (physician, nurse or clinical staff)
 - Fast measurement (a few minutes only)
 - Immediate results
 - Certified training & consumer service
-
- Clinically validated: 2,000+ peer-reviewed publications
 - Guidelines recommendations (EASL, AASLD, APASL, WHO)
 - Patented technology
 - Proprietary algorithms
-
- Adapted to most settings
 - Connectivity
 - Value-based
 - No consumable

OPTIONS

- Roll stand
- Dicom compatibility
- Maintenance contracts

RELATED PRODUCTS

- FibroMeter™ VCTE™
- FAST™ Score
- FibroScan® Gateway
- myFibroScan app

TECHNICAL PARAMETERS

- 📏 **Size:**
(H x L x D) 46 cm x 36 cm x 25 cm
- 📦 **Weight:**
10 kg approx. (including accessories)
- ⚡ **Power supply:**
HW version step 1: 100-240 V ~ 1.8 A 50/60 Hz
HW version step 2: 100-240 V ~ 1.2 A 50/60 Hz
- ➔ **Connection:**
Ethernet connector RJ45, 4 USB 2.0 ports, 2 probe connectors, 1 auxiliary output (footswitch compatible), DICOM compliant (optional)
- 🖱️ **Touch Screen:**
15-inch