


Edwards PASCAL Precision
Transcatheter Valve Repair System

PASCAL Precision System
Now Approved

Accurate,
intuitive control,
in your hands*



*Performance data on file and marketing evaluation.



Edwards

Designed for precise placement with accurate, intuitive control*

Advanced catheter and handle design facilitates smooth navigation and implant deployment*

Catheter control

Independent catheters allow for maneuvering within 3 planes

1

Secure stabilizers


Incremental movements and fine adjustments to advance, retract, and maneuver catheters*

2

Responsive catheter design

Optimized torque transfer to facilitate implant placement†

Laser-cut hypotube technology provides stability and balanced catheter flexibility



Hydrophilic coating and inter-catheter clearance for increased responsiveness



3


Ergonomic manipulation

Redesigned clasp sliders allow for both independent and simultaneous claspings, while implant actuation knob provides full implant configuration control‡

4

Simple release

Implant release mechanisms stay covered until the implant is ready to be deployed



Images are not actual size.

*Performance data on file and marketing evaluation.

†Design and performance data on file and marketing evaluation.

‡Design data on file and marketing evaluation.

Leading the way in MR and TR innovation

Edwards is committed to advancing the treatment of mitral regurgitation (MR) and tricuspid regurgitation (TR) for the patients who can be treated and to significantly improve their quality of life^{1,2}

2019



PASCAL Implant System for MR

Introduction of PASCAL repair system uniquely designed to enhance leaflet capture and close the coaptation gap¹

2020



PASCAL Implant System for TR

PASCAL Implant System for TR expanded indication for even challenging tricuspid anatomies*



PASCAL Ace Implant System for MR and TR

A narrower profile designed for ease-of-navigation in complex anatomies


2021



PASCAL Stabilizer Rail System

Enhanced stability to advance, retract, and torque catheters

2022



PASCAL Precision System

Delivery system catheters, handles, and stabilizers optimally designed for versatility, accuracy, and control

*Instructions For Use.



Innovation in your hands. [Learn more at Edwards.com/PASCAL](https://www.edwards.com/PASCAL)

Valve repair with the PASCAL implant system

Demonstrated consistent clinical results



97% freedom from recurrent MR at 2 years¹

78% of patients achieved MR≤1+

Study Design: Multicenter, prospective, single-arm study. Total patients: 124; Mean age: 74.9 years; Baseline MR Severity ≥3+: 100%.



100% of patients achieved 1 TR grade reduction at 1 year²

Study Design: Multicenter, prospective, single-arm study. Total patients: 65; Mean age: 78; NYHA functional class III or IV: 71%

Distinct implant designs, same functionality



PASCAL Ace Implant



PASCAL Implant



Access more details at [Edwards.com/PASCAL](https://www.edwards.com/PASCAL)

MR, mitral regurgitation; TR, tricuspid regurgitation; NYHA, New York Heart Association.

References:

1. Szerlip M, Spargias KS, Makkar R, et al. 2-Year outcomes for transcatheter repair in patients with mitral regurgitation from the CLASP study. *JACC Cardiovasc Interv.* 2021;14(14):1538-1548. 2. Hahn R. Transcatheter tricuspid valve repair: Clasp TR study one-year results. Presented at: LBT EuroPCR 2022; May 17-20, 2022; Paris, France. <https://course.pconline.com/en/vod/late-breaking-mitral-and-tricuspid-studies-and-registries>.

Medical device for professional use. For a listing of indications, contraindications, precautions, warnings, and potential adverse events, please refer to the Instructions for Use (consult [eifu.edwards.com](https://www.edwards.com) where applicable).

Edwards, Edwards Lifesciences, the stylized E logo, CLASP, Edwards PASCAL, PASCAL, PASCAL Ace, and PASCAL Precision are trademarks or service marks of Edwards Lifesciences Corporation. All other trademarks are the property of their respective owners.

© 2022 Edwards Lifesciences Corporation. All rights reserved. PP-EU-3759 v1.0

Edwards Lifesciences • Route de l'Etraz 70, 1260 Nyon, Switzerland • [edwards.com](https://www.edwards.com)



Edwards